CE 0483

N02969v2
CE0200 | CE0201 | CE2500

SERAFIT[®], SERAPID[®], SERAFIT[®] PROTECT

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POLYGLYCOLIC ACID SUTURE MATERIAL

- SERAFIT®
- SERAPID®
- SERAFIT® PROTECT

WARNINGS AND PRECAUTIONS

- Read the instructions for use in their entirety before using the product.
- The product should only be used by qualified personnel that is responsible to be familiar with both appropriate surgical techniques and anatomy prior to use.
- The product is delivered sterile. Sterilization method: ethylene oxide (SERAFIT®, SERAFIT® PROTECT), gamma irradiation (SERAPID®)

SERAFIT®
SERAFIT®
PROTECT

SERAPID®

- The product is intended for single use only.

- Do not re-sterilize or re-use the product as critical changes in performance and safety are possible.

- Do not use after expiration date.

- Do not use if the product shows signs of damage (e.g. cuts or kinks) or if the outer packaging or sterile packaging has been damaged or unintentionally opened prior to use.

- Do not use if the product is stored outside the recommended storage conditions.
- Use the product only for the intended purpose.

DESCRIPTION

SERAFIT[®], SERAPID[®] and SERAFIT[®] PROTECT are synthetic, absorbable, multifilament braided surgical suture materials. They are made from the homopolymer of glycolic acid and are coated. They are supplied dyed violet or undyed. More information is contained in the following table:

Product	Applied material	Thread structure	Dyestuff	Coating
SERAFIT [®]	Polyglycolic acid	Multifilament, braided	D+C violet No. 2 (Colour Index No. 60725) Content ≤ 0.2 wt % Also available undyed.	<ul style="list-style-type: none">• Calcium Stearate and Polycaprolactone or• Calcium Stearate and Poly-(glycolide-co-lactide) (30:70) Content ≤ 5 wt %
SERAPID [®]	Polyglycolic acid	Multifilament, braided	Not applicable (available only undyed).	<ul style="list-style-type: none">• Calcium Stearate and Polycaprolactone or• Calcium Stearate and Poly-(glycolide-co-lactide) (30:70) Content ≤ 5 wt %
SERAFIT [®] PROTECT	Polyglycolic acid	Multifilament, braided	D+C violet No. 2 (Colour Index No. 60725) Content ≤ 0.2 wt %	<ul style="list-style-type: none">• Calcium Stearate and Poly-(glycolide-co-lactide) (30:70) Content ≤ 5 wt %• Chlorhexidine diacetate• Content ≤ 60 µg/m

The suture material is available in various thread diameters (gauge sizes) and thread lengths. The thread length varies between 3 cm and 400 cm. The thread diameter varies between EP 0.4 and 8, respectively USP 8/0 and 6. The suture material is supplied without needle or as needle-thread combination with atraumatic needles. The needles are made of surgical steel and are coated with silicone. They are available in different sizes, shapes (e.g. circle, straight, semi-curved) and types (e.g. round bodied or reverse cutting).

For some articles of SERAFIT[®] the thread is provided in a suture guiding tube which is made of plastics (polyamide (PA), polybutylene terephthalate (PBT)).

Needles and suture guiding tubes are not intended to remain in the body.

The implantable suture material complies with the current edition of the European Pharmacopoeia (EP) monograph "Sutures, sterile synthetic absorbable braided (Fila resorbilia synthetica torta sterilia)".

INTENDED PURPOSE

The suture material is intended for sutures where synthetic absorbable multifilament thread material is required. The suture material is classified as an implantable device which is used invasively in the context of a surgical operation. Due to absorption it serves to approximate tissue temporarily during healing period and it loses tensile strength by hydrolysis.

INDICATION

Surgical sutures are indicated to attach tissue to tissue, tissue to foreign material and for ligatures which are a special type of suture used to close hollow organs.

In the choice and use of the device, due account has to be taken of the condition of the patient, the experience of the qualified personnel, the surgical technique, and the type and size of the wound.

SERAFIT[®] is used to attach tissue to tissue or tissue to foreign material. It is also used for ligation.

SERAPID[®] is used when a rapid absorption of the suture thread is required. It is used to attach tissue to tissue or tissue to foreign material.

SERAFIT® PROTECT is used when a preservation of the suture thread is required in order to prevent microbial colonization. It is used to attach tissue to tissue or tissue to foreign material. It is also used for ligation.

MODE OF ACTION

The performance of absorbable sutures is based on a mechanical mode of action. When introduced into a living organism, absorbable suture is absorbed by that organism. Due to hydrolysis, the suture material loses tensile strength over time.

SERAFIT®, SERAFIT® PROTECT: At least 50% of the initial tensile strength remains after 14 days (standardised in-vitro tests). Absorption is expected to be completed after about 90 days.

SERAPID®: At least 50% of the initial tensile strength remains after 5 days (standardised in-vitro tests). Absorption is expected to be completed after about 42 days.

SERAFIT® PROTECT: Chlorhexidine contained in the coating is a preservative agent to protect the suture material from microbial contamination.

APPLICATION NOTE

Standard knot-tying techniques should be used to ensure knot-holding security.

In order to prevent needle damage during suturing it is recommended that to grasp the needle at an area between the end of its first third (as seen from its end) and its midpoint.

In cases where the loop of articles provided in a suture guiding tube is not appropriately opened, it is recommended to push the tube slightly forwards.

For information about removal of the thread from the packaging, see corresponding chapter "instructions for removal from the packaging".

CONTRAINDICATIONS

Suture material made of polyglycolic acid is contraindicated for applications for which more prolonged or continuing stability of the suture material is required. The device is contraindicated for use in patients with allergies to any of the product components. The device is contraindicated for use in direct contact with the heart, the central circulatory system and the central nervous system.

SERAPID® is further contraindicated for use in neurosurgery or for ligatures.

SERAFIT® PROTECT is further contraindicated for use in microsurgery, orthopedics, ophthalmology and otology or at other sensitive organs and in newborns.

RESTRICTIONS AND LIMITATIONS OF USE / NOTES / INTERACTIONS

The requirements of the suture material and the techniques needed vary with the indication. In the choice of suture material, account should be taken of the characteristics of the suture material concerned.

The use of the device should be carefully considered in patients with delayed wound healing. Considerations should be taken regarding the use of absorbable sutures in tissues with poor blood supply as delayed absorption and suture extrusion may occur. The use of absorbable sutures may be inappropriate in elderly, malnourished or debilitated patients. Knots in interrupted suturing may impact the absorption time.

Like all foreign bodies, the device can give rise to calculus formation when in prolonged contact with salt-containing solutions (e.g. in the biliary or urinary tract).

Additional use of non-absorbable suture materials or suitable immobilization techniques may be required with sutures that are under tension or require greater support.

Use of suitable ancillary measures is required with sutures in infected tissue, as the multifilament structure of the thread can prolong infection. Because of the increased risk of infection associated with the multifilament structure of the thread, skin sutures should if possible be intracutaneous and the suture material, in particular the knots, should be placed

as deep as possible.

Skin sutures which remain in place longer than 7 days may cause local irritation and should be snipped off or removed.

Surgical instruments should be handled with care in order to avoid damage to the thread. Bending of needles or grasping them outside of the recommended area can result in needle breakage. Because of the risk of infection from needle pricks, needles should be handled with particular care.

UNDESIRABLE SIDE-EFFECTS

Dehiscence, infections, local irritation, wound inflammatory reactions, suture fistulas, suture erosion or wound healing disorder may occur.

Users and/or patients are encouraged to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)

After the launch of the European Database on Medical Devices (EUDAMED), the Summary of Safety and Clinical Performance (SSCP) of the medical device is located at the following web site: <https://ec.europa.eu/tools/eudamed>. Before that, the SSCP is made available to the public upon request.

STORAGE CONDITIONS

The storage conditions apply as indicated on the label.

DISPOSAL INFORMATION

Disposal of unused products or waste material with particular care due to the potential risk of contamination and/or injury. In all other cases observe the specific national regulations.

MEANING OF THE SYMBOLS

Symbols are explained at the end of the instructions for use.

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SERAG-WIESSNER GmbH & Co. KG