

**NASOPHARYNGEALTUBUS
AUS PVC**

- transparent
- steril



125410	05030997	12	2,5 mm	3,97 mm	60 mm	10
	05031034	14	3,0 mm	4,64 mm	72 mm	
	05031040	16	3,5 mm	5,30 mm	85 mm	
	05031057	18	4,0 mm	5,97 mm	105 mm	
	05031063	19	4,5 mm	6,29 mm	110 mm	
	05031086	20	5,0 mm	6,63 mm	115 mm	
	05031100	22	5,5 mm	7,29 mm	125 mm	
	05031117	24	6,0 mm	7,96 mm	130 mm	
	05031152	26	6,5 mm	8,63 mm	140 mm	
	05031229	28	7,0 mm	9,29 mm	155 mm	
	05031235	30	7,5 mm	9,96 mm	165 mm	
	05031241	32	8,0 mm	10,63 mm	170 mm	
	05031270	34	8,5 mm	11,29 mm	175 mm	
	05031301	36	9,0 mm	11,95 mm	180 mm	